

---

## Protein A Sepharose, From Staphylococcus aureus For Molecular Biology

**MB112**

---

### Product Identifier

HS Code : 3913 90 90  
Storage : 2 - 8°C  
Shelf life : 3 years

---

### Technical Specification

Appearance : White to off-white powder or tiny beads  
Solubility : Swells in water  
DNases : None detected  
RNases : None detected  
Loss on drying (at 105°C, 2 hr) : ≤ 4.0%

---

### Safety Information

UN No. : Not dangerous goods  
Class : -  
Packing Group : -  
RTECS : Not available  
WGK : 1