

---

## **p-Nitrophenyl phosphate disodium salt hexahydrate For Molecular Biology**

**MB124**

### **Product Identifier**

---

CAS No.	:	333338-18-4
EC No.	:	224-246-5
Molecular Formula	:	C <sub>6</sub> H <sub>4</sub> NNa <sub>2</sub> O <sub>6</sub> P.6H <sub>2</sub> O
Molecular Weight	:	371.14
HS Code	:	2919 90 90
Storage	:	2 - 8°C
Shelf life	:	4 years

### **Technical Specification**

---

Appearance	:	White to pale yellow crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Water (K.F.)	:	24.00 - 34.00%
Assay (NT)	:	min. 98.00%

### **Safety Information**

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1