
Starch hydrolysed

MB127

For Molecular Biology

Product Identifier

CAS No.	:	9005-25-8
EC No.	:	232-679-6
Molecular Formula	:	(C ₆ H ₁₀ O ₅) _n
Molecular Weight	:	(162.14) _n
HS Code	:	3505 10 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white powder or solid
Solubility	:	33.3 mg soluble in 1 mL of hot water
DNases	:	None detected
RNases	:	None detected
Sulphated ash	:	<= 1.00%
Loss on drying (at 105°C, 2 hr)	:	<= 20.0%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1