
Streptavidin, from *Streptomyces avidinii* For Molecular Biology

MB129

Product Identifier

CAS No.	:	9013-20-1
Molecular Weight	:	~60,000
HS Code	:	3504 00 99
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to off-white powder
Solubility	:	10 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Activity (Fibrinogen)	:	~14 units/mg protein

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1