

---

## Cellulose powdered, For Molecular Biology

**MB132**

### Product Identifier

---

CAS No.	:	9004-34-6
EC No.	:	232-674-9
HS Code	:	3912 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White to off-white powder
Solubility	:	Insoluble in water
pH (10% aqueous slurry at 25°C)	:	5.0 - 7.5
Identification test	:	Passes
Loss on drying (at 105°C, 2 hr)	:	<= 7.0%
Chloride (Cl)	:	<= 0.05%
Ash	:	<= 0.3%
Water soluble substances	:	<= 1.5%
Assay (FAS Titration, on dry basis)	:	97.00 - 102.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	FJ5691460
WGK	:	1