
Carboxymethylcellulose sodium salt, Medium viscosity (CMC), For Molecular Biology

MB138

Product Identifier

CAS No.	:	9004-32-4
HS Code	:	3912 31 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to light tan crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water forming colloidal solution
DNases	:	None detected
RNases	:	None detected
Viscosity (2% in water at 25°C)	:	400 - 800 cp

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1