
Xylan from beechwood, For Molecular Biology

MB141

Product Identifier

CAS No.	:	9014-63-5
EC No.	:	232-760-6
HS Code	:	2940 00 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to brown crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of 1M NaOH, may be turbid
DNases	:	None detected
RNases	:	None detected
Specific rotation	:	-100.00° to -50.00° (c = 1 in 5% NaOH at 20°C)
Loss on drying (at 105°C, 1 hr)	:	<= 10.00%

Safety Information

UN No.	:	Not Dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3