
Nonidet® P-40 (NP 40) substitute (Mixture of 15 homologues), MB143 For Molecular Biology


Product Identifier

CAS No.	:	9016-45-9
EC No.	:	500-024-6
Synonym	:	4-Nonylphenyl-polyethylene glycol; NP-40
HS Code	:	3402 42 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless viscous liquid
Solubility	:	33.3 mg soluble in 1 mL of water
pH (1% in water at 25°C)	:	5.00 - 8.00
DNases	:	None detected
RNases	:	None detected

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H411
Precautionary Statement(s)	:	P273- P280- P302 + P352- P305+P351+P338- P332 + P313
UN No.	:	3082
Class	:	9
Packing Group	:	III
RTECS	:	AX0247000
WGK	:	3