
Humic acid sodium salt

MB152

For Molecular Biology

Product Identifier

CAS No.	:	68131-04-4
EC No.	:	268-608-0
HS Code	:	3913 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Black crystals or powder or granules
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Loss on drying (at 105°C, 4 hr)	:	15.00 - 25.00 %
Assay (on carbon basis, as humic acid)	:	45.00 - 65.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	MT6550000
WGK	:	3