
Chelex[®]-100 Sodium form (50-100 Mesh), Cation exchange resin, For Molecular Biology

MB160

Product Identifier

CAS No.	:	11139-85-8
Molecular Weight	:	~3500
HS Code	:	3914 00 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to white crystals or beads
Solubility	:	Insoluble in water
Particle size	:	50 - 100 mesh (dry)
DNases	:	None detected
RNases	:	None detected
Cross linking	:	1% cross linked
Matrix	:	Styrene divinylbenzene
Wet bead size	:	300 - 1180 µm
Min. wet capacity	:	0.04 meq/mL
Operating pH range	:	4 - 14

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1