
Mineral oil, For Molecular Biology

MB161

Product Identifier

CAS No.	:	8042-47-5
EC No.	:	232-455-8
HS Code	:	2710 19 90
Storage	:	Below 30°C
Shelf life	:	2 years

Technical Specification

Appearance	:	Clear colourless liquid
Solubility	:	Insoluble in water
DNases	:	None detected
RNases	:	None detected
Density (at 25°C)	:	0.8390 - 0.8590 g/mL
Refractive index (n _{20/D})	:	~ 1.467

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	PY8047000
WGK	:	1