
Glycogen, From Oysters, Type II, For Molecular Biology

MB164

Product Identifier

CAS No.	:	9005-79-2
EC No.	:	232-683-8
Molecular Formula	:	(C ₆ H ₁₀ O ₅) _n
HS Code	:	3913 90 90
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to light beige powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water with opalescence
DNases	:	None detected
RNases	:	None detected
Total sugars (UV)	:	min. 90.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	MC2700000
WGK	:	1