

---

## Ficoll® 400, For Molecular Biology

**MB167**

---

### Product Identifier

CAS No.	:	26873-85-8
Molecular Weight	:	~ 400000
HS Code	:	3907 20 90
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to off-white solid or powder
Solubility	:	100 mg soluble in 1 mL of water to give clear to slightly hazy solution
DNases	:	None detected
RNases	:	None detected
Specific rotation	:	+50.00° to +59.00° (c = 1% in water at 20°C)

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1