

---

## Gelatin, Type A, For Molecular Biology

**MB169**

---

### Product Identifier

CAS No.	:	9000-70-8
EC No.	:	232-554-6
HS Code	:	3503 00 90
Storage	:	Below 30°C
Shelf life	:	3 years

---

### Technical Specification

Appearance	:	White to amber powder or granules or flakes
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Loss on drying (at 105°C, 2 hr)	:	≤ 16.0%
Total nitrogen	:	≥ 15.50%
Protein	:	70.00 - 90.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	LX8580000
WGK	:	1