
Agarose special, High EEO, For Molecular Biology (Nuclease and Protease free)

MB176

Product Identifier

| | | |
|------------|---|------------|
| CAS No. | : | 9012-36-6 |
| EC No. | : | 232-731-8 |
| HS Code | : | 3913 90 90 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|----------------------------|---|--|
| Appearance | : | White to off-white crystals or powder |
| Solubility | : | 33.3 mg soluble in 1 mL of boiling water |
| DNases | : | None detected |
| RNases | : | None detected |
| Protease | : | None detected |
| Gelling temperature | : | 34.5 - 37.5°C |
| Sulfate (SO ₄) | : | <= 0.2% |
| Electroendosmosis (-Mr) | : | 0.23 - 0.26 |
| Gelling strength (1% gel) | : | >= 750 g/cm ² |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | Not available |
| WGK | : | 1 |