

---

## TRIS acetate salt, For Molecular Biology

**MB190**

### Product Identifier

---

CAS No.	:	6850-28-8
EC No.	:	229-939-6
Molecular Formula	:	C <sub>6</sub> H <sub>15</sub> NO <sub>5</sub>
Molecular Weight	:	181.19
Synonym	:	Tris(Hydroxymethyl)aminomethane acetate salt
HS Code	:	2922 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White crystals or powder
Solubility	:	1000 mg soluble in 1 mL of water
pH (0.1 M in water at 25°C)	:	6.00 - 7.00
DNases	:	None detected
RNases	:	None detected
Assay (NaOH Titration, on dry basis)	:	99.00 - 102.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1