
ADP disodium salt, For Molecular Biology

MB191

Product Identifier

CAS No.	:	16178-48-6
EC No.	:	240-314-7
Molecular Formula	:	C ₁₀ H ₁₃ N ₅ Na ₂ O ₁₀ P ₂
Molecular Weight	:	471.16
HS Code	:	2934 99 90
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder.
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Water (K.F.)	:	<= 10.0%
Assay (Total Nitrogen/HPLC)	:	min. 95.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1