

Ammonium phosphate dibasic, For Molecular Biology

MB192

Product Identifier

CAS No.	:	7783-28-0
EC No.	:	231-987-8
Molecular Formula	:	(NH ₄) ₂ HPO ₄
Molecular Weight	:	132.06
HS Code	:	3105 30 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to white crystals or powder or chunks or solid
Solubility	:	100 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	7.80 - 8.50
DNases	:	None detected
RNases	:	None detected
Assay (HCl Titration)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1