

Light green SF yellowish

MB201

For Molecular Biology

Product Identifier

CAS No.	:	5141-20-8
EC No.	:	225-906-5
C.I. No.	:	42095
Molecular Formula	:	C ₃₇ H ₃₄ N ₂ O ₉ S ₃ Na ₂
Molecular Weight	:	792.85
HS Code	:	3204 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Red to dark brown or reddish purple or green crystals or powder or solid
Solubility	:	20 ppm in water yields clear solution
DNases	:	None detected
RNases	:	None detected
Absorption maxima	:	629 - 634 nm (20 ppm in water)
Loss on drying (at 110°C, 2 hr)	:	<= 12.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	BQ4900000
WGK	:	2