

---

## Sodium thiosulphate anhydrous

**MB210**

### For Molecular Biology

#### Product Identifier

---

CAS No.	:	7772-98-7
EC No.	:	231-867-5
Molecular Formula	:	Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub>
Molecular Weight	:	158.11
HS Code	:	2832 30 10
Storage	:	Below 30°C
Shelf life	:	4 years

#### Technical Specification

---

Appearance	:	Colourless or white crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Loss on drying (at 105°C, 2hr)	:	<= 2.00%
Assay (Iodimetry, on dry basis)	:	98.00 - 102.00%

#### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	XN6476000
WGK	:	1