
MOPSO, For Molecular Biology

MB214


Product Identifier

CAS No.	:	68399-77-9
EC No.	:	269-989-6
Molecular Formula	:	C ₇ H ₁₅ NO ₅ S
Molecular Weight	:	225.26
HS Code	:	2934 99 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
FTIR	:	Matches with the standard pattern
Melting range	:	275 - 280°C
Assay (NaOH Titration)	:	99.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1