

## Sodium cacodylate trihydrate, For Molecular Biology

**MB219**


### Product Identifier

CAS No.	:	6131-99-3
EC No.	:	204-708-2
Molecular Formula	:	C <sub>2</sub> H <sub>6</sub> AsO <sub>2</sub> Na.3H <sub>2</sub> O
Molecular Weight	:	214.03
HS Code	:	2931 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White crystals or powder or solid
Solubility	:	1000 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Assay (NT)	:	98.00 - 102.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301+H331- H350- H410
Precautionary Statement(s)	:	P201- P261- P273- P280- P301 + P310- P304 + P340- P501
UN No.	:	1688
Class	:	6.1
Packing Group	:	II
RTECS	:	CH7890000
WGK	:	3