
Dextran For Molecular Biology

MB221

Product Identifier

CAS No.	:	9004-54-0
EC No.	:	232-677-5
Molecular Formula	:	(C ₆ H ₁₀ O ₅) _n
Molecular Weight	:	~10,000
HS Code	:	3913 90 30
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH (1% in water at 25°C)	:	4.50 - 7.00
DNases	:	None detected
RNases	:	None detected
Average molecular weight	:	~ 10000

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	HH9230000
WGK	:	2