

Glutaraldehyde 25% w/w, For Molecular Biology

MB222


Product Identifier

CAS No.	:	111-30-8
EC No.	:	203-856-5
Molecular Formula	:	C ₅ H ₈ O ₂
Molecular Weight	:	100.12
HS Code	:	2912 19 90
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to yellow solution
Solubility	:	1 mL miscible in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Concentration (NaOH Titration)	:	24.50 - 30.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301- H314- H317- H330- H334- H410
Precautionary Statement(s)	:	P260- P280- P301 + P310- P303+P361+P353- P304 + P340- P305+P351+P338
UN No.	:	2927
Class	:	6.1(8)
Packing Group	:	II
RTECS	:	Not available
WGK	:	3