

Orcinol monohydrate, For Molecular Biology

MB242

Product Identifier

CAS No.	:	6153-39-5
EC No.	:	207-984-2
Molecular Formula	:	C ₇ H ₈ O ₂ .H ₂ O
Molecular Weight	:	142.15
HS Code	:	2907 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to tan to brownish-red crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of alcohol
DNases	:	None detected
RNases	:	None detected
FTIR	:	Matches with the standard pattern
Melting range	:	56 - 61°C
Assay (GC)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1