
Rubidium chloride, For Molecular Biology

MB243

Product Identifier

| | | |
|-------------------|---|------------|
| CAS No. | : | 7791-11-9 |
| EC No. | : | 232-240-9 |
| Molecular Formula | : | RbCl |
| Molecular Weight | : | 120.92 |
| HS Code | : | 2827 39 90 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|---------------------------|---|---|
| Appearance | : | White crystals or powder or chunks or solid |
| Solubility | : | 100 mg soluble in 1 mL of water |
| DNases | : | None detected |
| RNases | : | None detected |
| Assay (AT back Titration) | : | min. 99.50% |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | VL8575000 |
| WGK | : | 1 |