
Ribonucleic acid, From Torula yeast, For Molecular Biology

MB244

Product Identifier

CAS No.	:	63231-63-0
EC No.	:	277-256-7
Synonym	:	RNA
HS Code	:	2934 99 90
Storage	:	On receipt store at -20°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to tan powder or solid
Solubility	:	10 mg soluble in 1 mL of 0.1M acetate buffer, pH 5.0
DNases	:	None detected
RNases	:	None detected
Loss on drying (at 105°C, 2 hr)	:	<= 8.00%
Phosphorous (P)	:	>= 8.00%
Nitrogen content	:	>= 14.00%
Assay (Total nitrogen, on dry basis)	:	>= 86.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	VJ2249500
WGK	:	1