
Ribonucleic acid sodium salt, For Molecular Biology

MB245

Product Identifier

CAS No. : 73049-77-1
HS Code : 2934 99 90
Storage : 2 - 8°C
Shelf life : 3 Years

Technical Specification

Appearance : White to tan powder or solid
Solubility : 33.3 mg soluble in 1 mL of cold water
DNases : None detected
RNases : None detected
Loss on drying (at 105°C, 2 hr) : ≤ 10.0%
Nitrogen content : ≤ 13.0%

Safety Information

UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : Not available
WGK : 1