
X-Glucoside, For Molecular Biology

MB246

Product Identifier

CAS No.	:	15548-60-4
EC No.	:	239-603-0
Molecular Formula	:	C ₁₄ H ₁₅ BrClNO ₆
Molecular Weight	:	408.63
Synonym	:	5-Bromo-4-chloro-3-indolyl-β-D-glucopyranoside
HS Code	:	2938 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to off-white crystals or powder.
Solubility	:	33.3 mg soluble in 1 mL of dimethylformamide
DNases	:	None detected
RNases	:	None detected
Melting range	:	250 - 260°C
Assay (HPLC)	:	min. 97.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1