

---

**Paraffin wax (noncaking), For Molecular Biology**  
**Packed in SS container**

**MB265**

---

**Product Identifier**

CAS No.	:	8002-74-2
EC No.	:	232-315-6
HS Code	:	2712 20 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

**Technical Specification**

Appearance	:	Colourless or white translucent mass
Solubility	:	33.3 mg soluble in 1 mL of chloroform
DNases	:	None detected
RNases	:	None detected
Congearing Point	:	58 - 60°C

---

**Safety Information**

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RV0350000
WGK	:	nwg