
**Paraffin wax (noncaking), For Molecular Biology
Block packed in plastic container**

MB265B

Product Identifier

CAS No.	:	8002-74-2
EC No.	:	232-315-6
HS Code	:	2712 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white translucent blocks
Solubility	:	33.3 mg soluble in 1 mL of chloroform
DNases	:	None detected
RNases	:	None detected
Congearing Point	:	58 - 60°C

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RV0350000
WGK	:	nwg