
NADP disodium salt, For Molecular Biology

MB269

Product Identifier

CAS No.	:	24292-60-2
EC No.	:	246-129-8
Molecular Formula	:	$C_{21}H_{26}N_7Na_2O_{17}P_3$
Molecular Weight	:	787.37
Synonym	:	β -Nicotinamide adenine dinucleotide phosphate disodium salt
HS Code	:	2934 99 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to off-white powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Assay (HPLC)	:	min. 90.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1