

---

## Molecular sieve 3Å, Rods, For Molecular Biology

**MB280**

### Product Identifier

---

CAS No.	:	308080-99-1
Molecular Formula	:	$\text{KnNa}_{12-n}[(\text{AlO}_2)_{12}(\text{SiO}_2)_{12}] \cdot x\text{H}_2\text{O}$
Synonym	:	Potassium alumino silicate, Rods
HS Code	:	2842 10 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Beige coloured rods
Solubility	:	Insoluble in water
DNases	:	None detected
RNases	:	None detected

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1