
Calcium phosphate dibasic anhydrous

MB297

For Molecular Biology

Product Identifier

CAS No.	:	7757-93-9
EC No.	:	231-826-1
Molecular Formula	:	CaHPO ₄
Molecular Weight	:	136.06
HS Code	:	2835 25 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of dilute hydrochloric acid
DNases	:	None detected
RNases	:	None detected
Loss on ignition (at 500°C)	:	6.50 - 8.50%
Assay (EDTA Titration)	:	98.00 - 105.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1