
Chromium (III) chloride hexahydrate, For Molecular Biology MB298


Product Identifier

| | | |
|-------------------|---|--------------------------------------|
| CAS No. | : | 10060-12-5 |
| EC No. | : | 233-038-3 |
| Molecular Formula | : | CrCl ₃ .6H ₂ O |
| Molecular Weight | : | 266.45 |
| HS Code | : | 2827 39 90 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|-------------------|---|--|
| Appearance | : | Green to very dark green crystals powder |
| Solubility | : | 100 mg soluble in 1 mL of water |
| DNases | : | None detected |
| RNases | : | None detected |
| Assay (Iodometry) | : | 98.00 - 105.00% |

Safety Information

| | | |
|----------------------------|---|---|
| Hazard Pictogram(s) | : |  |
| Signal Word | : | Danger |
| Hazard Statement(s) | : | H290- H302- H314- H317- H411 |
| Precautionary Statement(s) | : | P260- P273- P280- P301+P330+P331- P303+P361+P353- P305+P351+P338- P310 |
| UN No. | : | 3260 |
| Class | : | 8 |
| Packing Group | : | III |
| RTECS | : | GB5450000 |
| WGK | : | 2 |