

## Paraffin liquid, light, For Molecular Biology

**MB352**


### Product Identifier

CAS No.	:	8012-95-1
EC No.	:	232-384-2
Synonym	:	Petrolatum liquid, light
HS Code	:	2710 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Colourless, transparent oily liquid
Solubility	:	33.3 mg soluble in 1 mL of ether
DNases	:	None detected
RNases	:	None detected
Density (at 25°C)	:	0.820 - 0.880 g/mL
Viscosity (at 30 ± 1°C)	:	10 - 40 cps (LV1 Spindle at 30 rpm)

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H304
Precautionary Statement(s)	:	P301 + P310- P331- P405- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	PY8030000
WGK	:	1