

---

## Span® 80, For Molecular Biology

**MB354**

### Product Identifier

---

CAS No.	:	1338-43-8
EC No.	:	215-665-4
Molecular Formula	:	C <sub>24</sub> H <sub>44</sub> O <sub>6</sub>
Molecular Weight	:	428.62
HS Code	:	3824 99 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Yellow to brownish yellow viscous liquid
Solubility	:	1 mL miscible in 1 mL of ethanol
DNases	:	None detected
RNases	:	None detected
Acid value (mg KOH/mg)	:	<= 8.00
Saponification value (mg KOH/mg)	:	145 - 160
Hydroxyl value (mg KOH/mg)	:	190 - 210

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	WG2932400
WGK	:	1