
L-Glutathione oxidized, For molecular Biology

MB357

Product Identifier

CAS No.	:	27025-41-8
EC No.	:	248-170-7
Molecular Formula	:	C ₂₀ H ₃₂ N ₆ O ₁₂ S ₂
Molecular Weight	:	612.63
HS Code	:	2930 90 99
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or solid or powder
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
Specific rotation	:	-106.0° to -92.0° (c = 2% in water at 20°C)
Water (K.F.)	:	<= 8.0%
Assay (NT/HPLC, on anhydrous basis)	:	98.00 - 101.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1