

2-Bromoethylamine hydrobromide, For Molecular Biology

MB384


Product Identifier

CAS No.	:	2576-47-8
EC No.	:	219-924-2
Molecular Formula	:	C ₂ H ₆ BrN.HBr
Molecular Weight	:	204.89
Synonym	:	2-Aminoethyl bromide hydrobromide
HS Code	:	2921 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to brown hygroscopic solid or crystals or powder or chunks
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
FTIR	:	Matches with the standard pattern
Melting range	:	169 - 178°C
Assay (AT)	:	98.00 - 106.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	KQ8225000
WGK	:	1