

MBLA003 : Thumb Press Dispensing Dropper/ Pasteur Pipette (Nuclease free), sterile

Product Number : Packing

MBLA003 : 1X50NO
: 1X500NO

Product Information

Product Code : MBLA003
Product Name : Thumb Press Dispensing Dropper/ Pasteur Pipette (Nuclease free), sterile

Technical Specification

Description : Thumb press dispensing dropper/ Pasteur pipettes are widely used in labs for sample and fluid transfer.

Clarity : semi-transparent

Color : Natural

Material : Polyethylene

Graduations : Yes (in 0.5ml increments)

Autoclavable : No

Capacity : 3.0 ml

Sterile : Yes

DNase/ RNase : None detected

Packaging : Individually packed

Use : Used to dispense liquid samples


Please refer disclaimer Overleaf.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



02/2026

PIMBLA003_0/0223

MBLA003-00

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg. office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com