

MBLA006 : Screw Cap Tubes with ‘O’ Ring

Product Number : Packing

MBLA006 : 1X1000NO

Product Information

Product Code : MBLA006

Product Name : Screw Cap Tubes with ‘O’ Ring

Technical Specification

Description : Self-standing, leak-proof, nuclease-free tube is a specialized laboratory consumable designed for secure, contamination-free handling and storage of sensitive biological samples.

Clarity : clear

Color : Natural

Material : Polyethylene

Graduations : No

Autoclavable : Yes

Capacity : 0.5 ml

Sterile : No

DNase/ RNase : None detected

Temperature Tolerance Limit : -80°C to 121°C

Use : Used to store liquid samples

Please refer disclaimer Overleaf.









Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.

Symbols:

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PIMBLA006

Rev No.: 01

Date of Issue.: 2025-12

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com