



# Technical Data

---

## MBLA010 Q4Pet Autoclavable Micropipette Capacity: 2 - 20 $\mu$ L

---

### Product Number Packing

MBLA010 : 1 No.

---

### Product Information

Product Code : MBLA010  
Product Name : Q4Pet Autoclavable Micropipette Capacity: 2 - 20  $\mu$ L

---

### Technical Specification

Description : Variable volume autoclavable microlitre pipettes wherein volumes can be continuously adjusted by rotating the plunger or thumb button. It has reverse pipetting technique. Each pipette has been tested and calibrated.  
Capacity :  
Use : 2 - 20  $\mu$ L  
: Used in laboratory to dispense liquid samples.

---

### Transport Information

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

#### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

---

Address : C-40, Road No.21/Y, MIDC, Wagle Ind. Area, Thane(W) - 400604, Maharashtra, India. Phone:+91-022-6147 1919, 6903 48 00 Fax:+91-022-6147 1920,6903 48 90 email : info@himedialabs.com. Visit us at our website: www.himedialabs.com