



Technical Data

MBLA012 Q4Pet Autoclavable Micropipette Capacity: 20 - 200 µL

Product Number Packing

MBLA012 : 1 NO

Product Information

Product Code : MBLA012

Product Name : Q4Pet Autoclavable Micropipette Capacity: 20 - 200 µL

Technical Specification

Description : Variable volume autoclavable microlitre pipettes wherein volumes can be continuously adjusted by rotating the plunger or thumb button. It has reverse pipetting technique. Each pipette has been tested and calibrated.

Capacity : 20 - 200 µL

Use : Used in laboratory to dispense liquid samples.

Transport Information

ADR/RID : Not Dangerous Goods

IMDG : Not Dangerous Goods

IATA : Not Dangerous Goods

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Address : C-40, Road No.21/Y, MIDC, Wagle Ind. Area, Thane(W) - 400604, Maharashtra, India. Phone:+91-022-6147 1919, 6903 48 00 Fax:+91-022-6147 1920,6903 48 90 email : info@himedialabs.com. Visit us at our website: www.himedialabs.com