

**MBLA017 : HiPer® Lock Microcentrifuge Tube, 2.0ml**

**Product Number : Packing**

MBLA017 : 1X500NO

**Product Information**

Product Code : MBLA017  
Product Name : HiPer® Lock Microcentrifuge Tube, 2.0ml

**Technical Specification**

Description : Micro centrifugal tubes with attached hinged cap. The cap has a double locking system which provides extra safety while handling precious/ hazardous samples. This unique capping system ensures safety of samples while processing. Frosted surface on the tube for easy labelling.

Clarity : Clear without any particles and fibers.

Color : Clear, natural

Material : Polypropylene

Graduations : Yes

Autoclavable : Yes

Height of the tube : Approx. 42mm

Outer Diameter : Approx. 10.75mm

Capacity : 2.0 mL

Temperature Tolerance

Range : -20°C to +121°C

Product alignment : Lid closure of the tube is tight.

Use : Suitable for Microcentrifugation

**Please refer disclaimer Overleaf.**

## Transport Information

---









ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

## Technical Assistance

---

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at [mb@himedialabs.com](mailto:mb@himedialabs.com).

### Symbols:

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PIMBLA017

Rev No.: 01

Date of Issue.: 2025-11

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

---

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.  
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com