

MBLA040 : Empty Tip box, 1ml

Product Number : **Packing**

MBLA040 : 1X5NO.

Product Information

Product Code : MBLA040

Product Name : Empty Tip box, 1ml

Technical Specification

Description : Autoclavable, stackable, micropipette tip box for 1000 µL tips, with 100 slots, hinged cover and lock system

Clarity : Clear without any particles and fibers.

Capacity to hold tips : 100 Units

Temperature tolerance test : Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.

Autoclavable : Yes

Dimensions

Length of box : 128mm

Width of box : 102mm

Total height of box : 97mm

Height of box without cap : 66mm

Length of inner tray : 125mm

Width of inner tray : 100mm

Height of inner tray : 10mm

Hole diameter : 8mm

Use : To store micropipette tips

Please refer disclaimer Overleaf.









Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.

Symbols:

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PIMBLA040

Rev No.: 00

Date of Issue.: 2025-02

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com