


# HiPer<sup>®</sup> Mini Vortexer MBLA047

Compact. Powerful. Affordable



## The Ideal Vortexer for Every Lab

The HiPer<sup>®</sup> Mini Vortexer combines powerful performance with a compact, ergonomic design making it the perfect choice for any laboratory environment. Whether you're mixing small volumes or working in tight spaces, this user-friendly vortex mixer ensures consistent results with a speed of 800-2,500 rpm.

## Key Features

- **Touch-Activated Mixing**  
Quick, responsive vortexing with minimal effort— just press and go.
- **800- 2,500 rpm Speed**  
Reliable, efficient mixing for everyday lab tasks.
- **Compact, Space-Saving Design**  
Perfect for crowded and tight workspaces.
- **Compatible with Tubes up to 15 mL**  
Versatile mixing for a wide range of common lab applications.
- **Dual Operation Modes**
  - **Touch Mode:** Activates when pressed
  - **Continuous Mode:** Runs non-stop for longer procedures

## Technical Specifications

Mode	Continuous/short-spin
Shaking movement	Orbital
Mixing orbit	3 mm
Speed	800~2500 rpm
Dimension (W×D×H) mm	74×91×84 mm
Input power	DC6 V
Weight	0.35 kg