

## Hi-FiBloE™ Nylon Membrane for Filtration

### MBM008- 50SH

HiMedia's Hi-FiBloE™ Nylon Membrane for Filtration can be used for biological applications like filtration of aqueous solutions and preparation of samples. These membranes are flexible, durable and tear resistant.

#### Applications:

Suitable for (biological analysis)

#### Specifications:


Colour: White

Diameter: 47mm

Pore size: 0.22µm

#### Storage:

HiMedia's Hi-FiBloE™ Nylon Membrane for Filtration can be stored at room temperature (15-25°C). Hi-FiBloE™ Nylon Membrane for Filtration is stable for 2 years when stored as directed.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,  
Reg. Off: Plot No. C-40, Road No.  
21Y, MIDC, Wagle Industrial Estate,  
Thane, (West) 400604, Maharashtra,  
INDIA. Web: [www.himedialabs.com](http://www.himedialabs.com)



02/2025

PIMBM008\_0/0222

MBM008-00

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Estate, Thane, (West) 400604, Maharashtra, INDIA. Customer Care No.: 00-91-22-6116 9797 Tel : 00-91-22-6147 1919, 6903 4800 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)



#### Registered Office

#### HiMedia Laboratories Pvt Ltd.

Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Estate,  
Thane, (West) 400604, Maharashtra, INDIA.  
Customer Care No.: 00-91-22-6116 9797  
Tel : 00-91-22-6147 1919, 6903 4800

Fax : 6147 1920  
Web : [www.himedialabs.com](http://www.himedialabs.com)  
Email : [info@himedialabs.com](mailto:info@himedialabs.com)  
[mb@himedialabs.com](mailto:mb@himedialabs.com)

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein