


MBT057

dTTP Solution, 100mM, 0.25ml

Product Name	Product Code	Kit Packing
dTTP solution, 100 mM, 0.25 ml	MBT057-0.25ML	0.25 ml



Formula: C₁₀H₁₄N₂Na₃O₁₄P₃

Description:

dTTP (2'-deoxythymidine 5'-triphosphate) is supplied as 100 mM aqueous solution at pH 7.0.

Features:

- Greater than 99% purity confirmed by HPLC
- Free of Human and *E.coli* DNA
- Highly stable

Applications:

For use in PCR, sequencing, nick translation, fill-in, cDNA synthesis, and TdT tailing reaction, as well as for dilution of radiolabeled dNTPs.

Storage conditions:

dTTP solution should be stored at -20°C. When stored under the recommended conditions, the product is stable for 2 years.

Quality Control:

Functionally tested in PCR


Greater than 99% purity confirmed by HPLC

Free of endo-, exo-deoxyribonuclease, ribonuclease and nicking activities

Technical Assistance

At HiMedia's, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.

Please refer disclaimer Overleaf.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
www.himedialabs.com



01/2025

PIMBT057_O/0122

M8T057-05

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com