



## Buffered HiCynth™ Peptone Water

MCD614

Buffered HiCynth™ Peptone Water is a pre-enrichment medium used for pre-enrichment of injured *Salmonella* species from foods prior to selective enrichment and isolation.

### Composition\*\*

Ingredients	Gms / Litre
HiCynth™ Peptone No.1*	10.000
Sodium chloride	5.000
Disodium phosphate	3.500
Monopotassium phosphate	1.500
Final pH ( at 25°C)	7.2±0.2

\*\*Formula adjusted, standardized to suit performance parameters

\*Chemically defined peptone

### Directions

Suspend 20.00 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense in 50 ml amounts. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

### Principle And Interpretation

Buffered HiCynth™ Peptone Water is a pre-enrichment medium designed to help recovery of sub-lethally damaged Salmonellae before transfer to a selective medium. It is prepared by completely replacing animal or vegetable peptones with chemically defined peptones to avoid BSE/TSE risks associated with animal peptones. This pre-enrichment medium is free from inhibitors and is well buffered and provides conditions for resuscitation of the cells that have been injured by processes of food preservation. It was noted by Edel and Kampelmacher (1) that sub-lethal injury to *Salmonella* may occur due to food preservation techniques involving heat, desiccation, high osmotic pressure, preservatives or pH changes. Buffered HiCynth™ Peptone Water during the pre-enrichment period helps in recovery of injured cells that may be sensitive to low pH (2). This is particularly important for vegetable specimens, which have low buffering capacity. This medium can be used for testing dry poultry feed (3). In a survey involving isolation of Salmonellae from meat that had been artificially contaminated with sub-lethally injured organisms, pre-enrichment in Buffered Peptone Water at 37°C for 18 hours before selection in Tetrathionate Brilliant Green Bile Broth showed superior results compared with direct selection method. Lactose Broth is frequently used as a pre-enrichment medium but it may be detrimental to recovery of Salmonellae (4).

The media contain HiCynth™ Peptone No.1 as a source of carbon, nitrogen, vitamins and minerals. Sodium chloride maintains the osmotic balance and phosphates buffer the medium. The broth is rich in nutrients and produces high resuscitation rates for sublethally injured bacteria and supports intense growth. The phosphate buffer system prevents bacterial damage due to changes in the pH of the medium.

Inoculate 10 grams specimen in 50 ml of these media and incubate at 35-37°C for 18 hours. Transfer 10 ml from this medium to 100 ml of Tetrathionate Broth and incubate at 43°C for 24 - 48 hours and then subculture on selective plating media.

Examine the plates for characteristic *Salmonella* colonies.

### Quality Control

#### Appearance

Cream to yellow homogeneous free flowing powder

#### Colour and Clarity of prepared medium

Light yellow coloured, clear solution without any precipitate

#### Reaction

Reaction of 2.0% w/v aqueous solution at 25°C. pH : 7.2±0.2

#### pH

7.00-7.40

#### Cultural Response

Please refer disclaimer Overleaf.

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.(Recovery is carried out using XLD Agar,M031).

### Cultural Response

Organism	Inoculum (CFU)	Growth	Recovery
<b>Cultural Response</b> <i>Salmonella Enteritidis</i> ATCC 50-100 13076		good-luxuriant	>=50%
<i>Salmonella Typhi</i> ATCC 6539	50-100	good-luxuriant	>=50%
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	good-luxuriant	>=50%

### Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

### Reference

1. Edel and Kampelmacher, 1973, Bull. W.H.O., 48:167.
2. Sadovski, 1977, J. Food Technol., 12:85.
3. Juven, Cox, Bailey, Thomson, Charles and Schutze, 1984, J. Food Prot., 47:299.
4. Angelotti, 1963, Academic Press, New York, N.Y.

Revision : 00 / 2015

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.