

### Crime Scene Investigation Kit

<u>Product Name</u>	<u>Product Code</u>	<u>Kit Packing</u>
Crime Scene Investigation Kit	ML201-1PK	1pk

**Introduction:** HiMedia's Crime Scene Investigation Kit has been designed exclusively to help forensic investigators at crime location.

#### Kit Contents

S.No	Product Code	Product Name	Quantity
			1PK
1	CO002	Triclogel (50ml)	1 No
2	C0100	SteriSwift™ Disinfectant wipes	1 pack
3	LA227	Disposable Surgeon's Kit (XL)	3 No
4	LA892	Latex Examination Gloves, Large Size	4 Pairs
5	LA706	Stainless Steel Forceps, blunt, 6 inch	2 No
6	LA708	Stainless Steel Forceps, blunt, 10 inch	2 No
7	LA709	Stainless Steel Forceps, pointed, 6 inch	2 No
8	LA711	Stainless Steel Forceps, pointed, 10 inch	2 No
9	DS0593	Scissor 6 inch, Sharp	1 No
10	DS0594	Scissor 8 inch, Sharp	1 No
11	DS0595	Scissor 6 inch, Blunt	1 No
12	DS0596	Scissor 8 inch, Blunt	1 No
13	PW1280	Sterile Cotton Swab	20 No
14	ML189	Blood Identifier Kit	2X25R
15	LA699	Cryogenic Permanent Marker, BLACK	1 No
16	LA700	Cryogenic Permanent Marker, BLUE	1 No
17	LA701	Cryogenic Permanent Marker, RED	1 No
18	DS0598	UV protectant Safety goggle	1 No
19	DS0599	Blue pen	2 No
20	LA769	Sterile Scalpel Blade	20 No
21	LA829	Blade holder	2 No
22	DS0602	Butter Paper	15 No
23	ML024	Molecular Biology Grade Water	30ml
24	DS0603	Filter Paper	10 No

Please refer disclaimer Overleaf.

### Storage conditions:


Crime Scene Investigation Kit has to be stored at 15-25 °C.

### Safety Information

The Crime Scene Investigation Kit is for laboratory use only, not for drug, household or other uses. Please refer the Material Safety Data Sheet (MSDS) for information regarding hazards and safe handling practices.

### Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at [mb@himedialabs.com](mailto:mb@himedialabs.com).



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,  
Reg. Off: Plot No. C-40, Road No. 21Y,  
MIDC, Wagle Industrial Area, Thane,  
(West) 400604, Maharashtra, INDIA.  
Web: [www.himedialabs.com](http://www.himedialabs.com)



02/2025

PIML201\_0/0222

ML201-01

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.  
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)